

A bloody companion.

Reach highest smear consistency with Etzel I.

- **Reproducible blood smears** through automatic control
- **Consistently high smear quality**
- **Continuously adjustable** smear speed
- **Variable length and thickness** of blood smears
- **Affordable Technology** - "Made in Germany"

Technical specifications

Slide Capacity	2 Slides
Power requirements	100-240 V AC
Output voltage	DC 12 V / 1 A
Dimensions (Width x Height x Depth), [cm (IN.)]	12,5 x 11,5 x 17 (5 x 4,5 x 6,5)
Weight, [kg (lb)]	1,5 (3,3)



Bio-Gram Diagnostics GmbH
Technologiezentrum Ludwigshafen
Donnersbergweg 1
D-67059 Ludwigshafen, Germany

sales@bio-gram.de