

A part of the German Nibelungen Legend.

Alberich I is a miniaturized but powerful bench top mixer which is recommended for laboratories, clinics, hospitals and pharmacies.

- **Force activation**
- **Minimalistic Design**
- **Suitable for different container types**
- **Continuously adjustable** rotation speed (3000 - 4500 rpm)
- **For liquid containers up to 50 ml**
- **Stops automatically** when vertical force is decreased
- **Suitable power adapter for your country** available

Technical Specifications

Power requirements	100-240 V AC
Output Voltage	DC 12 V / 1 A
Dimensions (Width x Height x Depth), [cm (IN.)]	10,5 x 7 x 10,5 (4 x 3 x 4)
Weight, [kg (lb)]	0,5, (1,1)
Recommended ambient temperature	5 °C bis 35 °C
Rotation speed	3000-4500 rpm
Max. loading weight	0,1 kg
Protection Class	IP53
Rotating diameter	4,5 mm



Bio-Gram Diagnostics GmbH
Technologiezentrum Ludwigshafen
Donnersbergweg 1
D-67059 Ludwigshafen, Germany

sales@bio-gram.de